



Manufacturing Equipment for Abatement and Restoration Industry Since 1988

PEROx 1200 Air Purifier Module for Medical Facilities

Highly Effective on MRSA and SARS virus; Destroys Microbe, Virus, Odor and VOCs

The PEROx air purifier module utilizes IHP (Ionized-Hydro-Peroxide) oxidizing processes to effectively destroy a wide range of harmful pathogens such as MRSA, SARS viruses, strong odors, VOCs and other bio-contaminants. Oxides are generated from moisture and oxygen in the environment and unused oxides convert back to O2 and H2O. PEROx air modules can safely and effectively sanitize the air and surrounding hard surfaces in critical care environments, such as hospitals, labs or nursing homes.

The PEROx air purifier module attaches directly to the exhaust end of the OmniAire 1200PAC portable HEPA air cleaners. Depending on size of the space treated, duration and required level of sanitation required, the PEROx offers LOW/MED/HIGH levels of IHP output. The LOW output is used in small and medium sized, occupied rooms and the HIGH output is used for sanitizing hospital rooms before and after occupancy

SPECIFICATIONS:

- Controls:** ON-OFF-TIMER Switch, Hour Meter, 0-12 hr timer; LOW/MED/HIGH OUTPUT selector switch
- Power requirements:** 115 VAC/60 Hz/4.5 Amps
- Purifier:** x2 PHI cell- Ionized Hydro-Peroxide generators
- Airflow:** 200-1100 CFM From 1200 PAC machine
- Pre-filter:** Pleated Filter, MERV-11 efficiency
- Housing:** Powder coated Aluminum, antique white
- Weight with Filters:** 24 lbs.
- Operation environment:** Non-hazardous environment, 104° F temperature max.
- Dimensions:** 8” high housing.

ORDERING INFOPART

OmniAire PEROx 1200PAC.....	OAHP1200
Primary/Secondary Filter	OPF1616
Quick Clamp, 10”.....	QCW12
Flexible Duct, 10”x25’	OAD10
5” Replacement Generator Cell.....	EHPG-2
9” Replacement Generator Cell.....	EHPG-3

Omnitec Design, Inc.
4640 Campus Place, Suite 110 Mukilteo, WA 98275
Visit us @ www.omnitecdesign.com
Office: 425-290-3922